

Integra®

DuraSeal® Dural Sealant System
Assembly & Use

Limit uncertainty by using the quick and easy to use DuraSeal® dural sealant system to ensure a watertight closure.



INTEGRA®
LIMIT UNCERTAINTY

Dual Liquid Applicator

Kit Components



Powder Vial



Blue Precursor Syringe



Clear Precursor Syringe



Plunger cap



Syringe holder








Spray Tips



Dual liquid applicator

Step 1 • Preparation



- 1-1** Depress threaded fitting of vial cap (Figure 1-1a). Ensure red line is no longer visible (Figure 1-1b).
- 1-2** Remove syringe cap from blue precursor syringe. Screw blue precursor syringe onto powder vial (Figure 1-2a), then inject syringe contents into vial (Figure 1-2b).
- 1-3** Gently shake vial/syringe assembly until powder is completely dissolved. The solution will turn blue (Figure 1-3).
- 1-4** Invert vial/syringe assembly, then draw vial contents back into syringe (Figure 1-4a). Unscrew syringe from vial and discard vial (Figure 1-4b).
- 1-5** Remove syringe cap from clear precursor syringe. Prior to applicator assembly, eliminate bubbles from both syringes and ensure precursor volume in two syringes is equal (Figure 1-5).






Step 2 • Assembly

- 2-1** Screw clear and blue precursor syringes onto applicator (Figure 2-1).
- 2-2** Slide syringe holder (Figure 2-2a) along syringe barrels until it fits snugly against syringe flanges. Attach plunger cap (Figure 2-2b) to syringe plungers.
- 2-3** Attach a spray tip to the Dual Liquid Applicator (Figure 2-3).



Final step • Application

- 3-1** Prepare application site by removing all blood clots and fluid.
- 3-2** Do not prime Dual Liquid Applicator prior to use, as plugging may result.
- 3-3** Position Dual Liquid Applicator tip (Figure 3-3) approximately 2 cm from target site. Rapidly depress syringes to ensure precursor mixing and spray-like application of gel.
- 3-4** Apply sealant to create 1-2 mm maximum thickness.




Extended Tip Applicator

Kit Components



Step 1 • Preparation



- 1-1** Depress threaded fitting of vial cap (Figure 1-1a). Ensure red line is no longer visible (Figure 1-1b).
- 1-2** Remove syringe cap from blue precursor syringe. Screw blue precursor syringe onto powder vial (Figure 1-2a), then inject syringe contents into vial (Figure 1-2b).
- 1-3** Gently shake vial/syringe assembly until powder is completely dissolved. The solution will turn blue (Figure 1-3).
- 1-4** Invert vial/syringe assembly, then draw vial contents back into syringe (Figure 1-4a). Unscrew syringe from vial and discard vial (Figure 1-4b).
- 1-5** Remove syringe cap from clear precursor syringe. Prior to applicator assembly, eliminate bubbles from both syringes and ensure precursor volume in two syringes is equal (Figure 1-5).






Step 2 • Assembly

- 2-1** Screw clear and blue precursor syringes onto applicator (Figure 2-1a).
- 2-2** Slide syringe holder (Figure 2-2 ①) along syringe barrels until it fits snugly against syringe flanges. Attach plunger cap (Figure 2-2 ②) to syringe plungers.
- 2-3** Adjust malleable applicator shaft to improve access or visualization (Figure 2-3).



Final step • Application

- 3-1** Prepare application site by removing all blood clots and fluid.
- 3-2** While in surgical field, whenever anatomically possible, briefly spray sealant on gauze and without interrupting flow move to the target site (Figure 3-2).
- 3-3** Position applicator tip 2-4 cm from the target site, depress syringes using strong, even pressure (Figures 3-3). Using a continuous motion, apply an even, approximately 1-2 mm thin coating of hydrogel.



MicroMyst® Applicator





Kit Components



MicroMyst® Applicator*
(to be ordered separately)

Step 1 • Preparation

- 1-1** Depress threaded fitting of vial cap (Figure 1-1a). Ensure red line is no longer visible (Figure 1-1b).
- 1-2** Remove syringe cap from blue precursor syringe. Screw blue precursor syringe onto powder vial (Figure 1-2a), then inject syringe contents into vial (Figure 1-2b).
- 1-3** Gently shake vial/syringe assembly until powder is completely dissolved. The solution will turn blue.
- 1-4** Invert vial/syringe assembly, then draw vial contents back into syringe (Figure 1-4a). Unscrew syringe from vial and discard vial (Figure 1-4b).
- 1-5** Remove syringe cap from clear precursor syringe. Prior to applicator assembly, eliminate bubbles from both syringes and ensure precursor volume in two syringes is equal (Figure 1-5).







*MicroMyst® applicator requires an open air source to operate - used in conjunction with Flow Regulator
IMPORTANT: Please refer to package insert for complete instructions, contraindications, warnings and precautions.




Step 2 • Assembly

- 2-1** Ensure the Flow Regulator is connected to a compressed air or Nitrogen source (Figure 2-1). Confirm the nitrogen source is set to 50 - 200psi (3.45 - 13.8 Bar).
- 2-2** Connect MicroMyst® applicator line to Flow Regulator (Figure 2-2a). Ensure air is flowing through MicroMyst® applicator (Figure 2-2b).
- 2-3** Position syringe holder over applicator fittings (Figure 2-3a). Screw clear and blue precursor syringes onto applicator (Figure 2-3b).
- 2-4** Slide syringe holder (Figure 2-4 ①) along syringe barrels until it fits snugly against syringe flanges. Attach plunger cap (Figure 2-4 ②) to syringe plungers.
- 2-5** Adjust malleable applicator shaft to improve access or visualization (Figure 2-5).



Final step • Application

- 3-1** Prepare application site by removing all blood clots and fluid.
- 3-2** While in surgical field, whenever anatomically possible, briefly spray sealant on gauze and without interrupting flow move to the target site (Figures 3-2).
- 3-3** Position applicator tip 1-4 cm from the target site, depress syringes using strong, even pressure (Figures 3-3). Using a continuous motion, apply an even, approximately 1-2 mm thin coating of hydrogel.



DuraSeal® Dural Sealant System: the PEG* Hydrogel Technology for a Watertight Dural Closure

| STRENGTH ² | BIOCOMPATIBILITY ^{1,3} | VISIBILITY | SPEED ^{1,4} | CONVENIENCE |
|-------------------------------------------------------------------------|---------------------------------------------------|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|-----------------------------------------|
| Tissue adherence and cohesive strength to withstand critical pressures. | Biocompatible PEG (polyethylene glycol) hydrogel. | Distinctive blue colorant provides visualization to assess sealant coverage and thickness. | Prepared in less than two minutes. When applied, forms a watertight seal in seconds. | Single kit, stored at room temperature. |

Indications:

The DuraSeal® dural sealant system is intended for use as an adjunct to standard methods of dural repair, such as sutures, to provide watertight closure.

The Extended Tip Applicator is intended for use in the simultaneous delivery of two non-homogenous solutions onto a surgical site.

The MicroMyst® Applicator is intended for use in the delivery of two non-homogenous solutions onto a surgical site.

The Flow Regulator is intended to provide pressurized gas (air or nitrogen) to gas-assisted applicators.

Contraindications:

Do not apply the DuraSeal® Dural Sealant in abdominopelvic surgical procedures for use as a sealant or adhesion barrier.

Do not use Extended Tip Applicator, MicroMyst® Applicator and Flow Regulator for other indications than ones provided in the instructions for use.

1. Cosgrove GR et al. "Safety and efficacy of a novel polyethylene glycol hydrogel sealant for watertight dural repair". J Neurosurg 106: 52-58, 2007
 2. In-vitro Product Comparison Study of Wound Healing Sealants. Report no: R090417B (Cyanta Report)
 3. Safety Testing for 4a20kSG-trilysine sealant with BHT. Report no: ER1105, page 2-9
 4. Delashaw JB. et al. "Reconstruction After Posterior Cranial Fossa Surgery: Application of a Synthetic Tissue Sealant to Augment Dural Closure." Case Report 2009
- * Polyethylene Glycol

Integra Contact - Sales & Marketing EUROPE, MIDDLE-EAST and AFRICA

Integra LifeSciences Services (France) SAS
Immeuble Séquoia 2 97 allée Alexandre Borodine
Parc technologique de la Porte des Alpes
69800 Saint Priest FRANCE
Phone +33 (0)4 37 47 59 00
Fax +33 (0)4 37 47 59 99
emea.info@integralife.com

LATIN AMERICA, ASIA PACIFIC and CANADA

Integra LifeSciences Corporation
311 Enterprise Drive
Plainsboro, New Jersey 08536 USA
Phone +1 609 936 5400 ■ +1 609 750 4259 (fax)
lapac_cs@integralife.com

Integra Contact - Questions & Information

Europe, Middle East and Africa: +33 (0)4 37 47 59 50 ■ +33 (0)4 37 47 59 25 (fax) ■ csemea@integralife.com
Benelux: +32 (0)2 257 4130 ■ +32 (0)2 253 2466 (fax) ■ custsvcbenelux@integralife.com
United Kingdom: +44 (0)1 264 345 780 ■ +44 (0)1 264 363 782 (fax) ■ custsvcs.uk@integralife.com
Switzerland: +41 (0)2 27 21 23 30 ■ +41 (0)2 27 21 23 99 (fax) ■ custsvcsuisse@integralife.com
Germany: +49 (0)2102 5535 6200 ■ +49 (0)2 102 5536 636 (fax) ■ custsvcgermany@integralife.com
Australia: +61 (0) 3 8540 0400 ■ +61 (0) 3 9540 0004 (fax) ■ orders@integra-LS.com.au
Latin America: +1 609 936 5400 ■ +1 609 750 4259 (fax) ■ lapac_cs@integralife.com
Asia Pacific: +1 609 936 5400 ■ +1 609 750 4259 (fax) ■ lapac_cs@integralife.com
New Zealand: +64 (0) 9 447 3460 ■ +64 (0) 9 447 3440 (fax) ■ customerservice@integra-LS.com.au
Canada: +1 800 263 6693 or +1 905 632 788 ■ +1 905 632 7938 (fax) ■ custserv.canada@integralife.com
integralife.eu - integralife.com

Contact - Product Ordering

integralife.eu/duraseal-order/



Covidien LLC
15 Hampshire Street
Mansfield, MA 02048 ■ USA



Covidien Ireland Limited
IDA Business and Technology Park
Tullamore ■ Ireland



Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Please read carefully the instructions for use.
- Non contractual document. Integra reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.

Additional information for EMEA Customers only:

Products mentioned in this document are CE class IIa and III devices. Please contact Integra customer service should you need any additional information on devices classification. All the medical devices mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED". Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation in the United States and/or other countries. DuraSeal and MicroMyst are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.

Document for use outside the USA only. ©2014 Integra LifeSciences Corporation. All rights reserved. Last modification: 03/2014. 0237673