


Integra CRW Stereotactic Surgical Products

Head Rings and Accessories

UCHRA



Universal Compact Head Ring Assembly

CT and MR scanner compatible when used with Integra Head Ring Screws; for use with BRWLF and LLO1

Includes

- Universal Head Ring
- Adapter Plate
- 2 Anterior and 2 Posterior Head Posts
- Cross Bar
- Cross Bar Screws
- Ear Bar Assembly
- 2 Head Ring Wrenches
- Cleaning Tap
- Sterilization Case
- User's Manual

UCHRAP



Arc Adapter Plate

The Adapter Plate allows attachment of the BRWLF and CRWPRECISE to the UCHRA

UCHRCASE



Storage and Sterilization Case

Only compatible with the UCHRA

Integra CRW Stereotactic Surgical Products

Head Rings and Accessories

UCHRP



Post and Cross Bar Kit

Includes

- 1 Cross Bar—Allows alternative positioning of head ring screws
- 2 Anterior/2 Posterior Posts
- 2 Cross Bar Screws

Available Separately

- UCHRCB1—Head Ring Cross Bar
- UCHRPA—1 Anterior Post (short)
- UCHRPP—1 Posterior Post (long)
- UCHRCBS—Head Ring Cross Bar Attachment Screw

UCHREBA




Universal Compact Head Ring Ear Bar Kit

Includes (2 each)

- Head Ring Supports
- Ear Posts

UCHREB - Ear Bars, Includes 1 each - ear post

HRW



Head Ring Wrench

A universal wrench used on many of the screws found in the Integra CRW system including the Head Ring drive's posts, and screws.

UCHRHK

Universal Compact Head Ring Hardware Kit

Includes

- Adapter Plate Attachment Screws (4)
- Intubation Hoop Attachment Screws (2)
- Head Post Attachment Screws (4)
- Head Post Wedges (4)
- Ear Bar Assembly Attachment Screws (2)
- Ear Bar Assembly Nylon Thumb Screws (2)

Integra CRW Stereotactic Surgical Products

Head Rings and Accessories

HRAIM



Intubation Head Ring Assembly

CT scanner compatible; for use with BRWLF

Includes

Intubation Head Ring

4 Head Ring Drives

4 Head Ring Posts


Head Ring Wrench

Cleaning Tap

Head Ring Positioner

2 Conical Sliding T-bolt Screws


HRD



Head Ring Drive

Head Ring Posts are mounted on Head Ring Drives. Head Ring Drives are inserted into HRAIM.

HRDXS



Head Ring Drive Extender/ Shortener

Lengthens or Shortens the working length of the Head Ring Posts to allow treatment of low tumors or a wide range of head sizes.

HRP



Head Ring Post

Only compatible with the HRAIM

TAP



Cleaning Tap for HRP

(Not for use with Aluminum Posts)

Integra CRW Stereotactic Surgical Products

Head Rings and Accessories

HRKTP



Head Ring Positioner

Includes (1 each)

Positioner Ring (not sold separately)

Velcro Strap (Also available in boxes of 5, use part no. GTCSV5)

CSS



Conical T-bolt Screw

Includes

CSS (1 short)–Used to attach the HRAIM and UCHRA Head Rings to the clamping plate and MAYFIELD® adapter

HREBA



Intubation Head Ring Ear Bar Kit

Includes (2 each)

Head Ring supports

Ear Posts

HRIMS



Head Ring Sector

Converts the assembly to the closed ring configuration

HRIMTS



T-handle Screw

For Intubation Hoop