

Thoracic Catheters, Trocars and Umbilical Catheters

Trocar Catheters

- Thermosensitive PVC with satin slip surface catheter with x-ray line and latex free

Code	Description	Diameter Fr	Length	Sales Unit
211-10	Suture rib trocar catheter	10	22cm	10
211-12	Suture rib trocar catheter	12	25cm	10
211-16	Suture rib trocar catheter	16	25cm	10
211-20	Suture rib trocar catheter	20	40cm	10
211-24	Suture rib trocar catheter	24	40cm	10
211-28	Suture rib trocar catheter	28	40cm	10
211-32	Suture rib trocar catheter	32	40cm	10



Argyle™ Umbilical Catheters

- Single-use, in transparent polyurethane, resistant to lipids, soft and thermosensitive, latex free
- Indicated for measuring pressure, taking samples, infusions and transfusions
- For continuous or intermittent access to the umbilical artery or vein
- Single lumen and dual lumen

Code	Description	X-Ray Opaque Lines	Dimensions Fr	Dimensions Length	Sales Unit
Single Lumen Catheters					
8888160325	Polyurethane	3	2.5	30.5cm	10
8888160333	Polyurethane	3	3.5	38cm	10
8888160341	Polyurethane	3	5.0	38cm	10



Code	Lumens	Inside Diameter G	Dimensions Fr	Dimensions Length	Sales Unit
Dual Lumen Catheters					
8888160531	2	20-23	3.5	38cm	5
8888160556	2	18-21	5	38cm	5
8888160648	3	18-21-21	5	38cm	5
8888160663	3	15-18-18	8	38cm	5



Cardiothoracic Accessories

Argyle™ Double Lumen Left Ventricular Vent Catheter

- Single-use catheter for left ventricular venting, in thin-walled, thermosensitive PVC, latex free

Code	Dimensions Fr	Dimensions Outside Diameter	Sales Unit
8888590208	14	4.7mm	10
8888590216	18	6mm	10
8888590224	20	6.7mm	10



Cardiothoracic Accessories

Argyle™ Aortic Perfusion Cannula



- PVC, single-use, thin-walled internal lumen, for maximum flow with minimum trauma, latex free

Code	Description	Dimension Fr	Tip Length	Sales Unit
8888591008	Angled cannula	10	19	10
8888591123	Angled cannula	12	20	10
8888591024	Angled cannula	14	22	10
8888591149	Angled cannula	16	23	10
8888591040	Angled cannula	18	25	10
8888591065	Angled cannula	21	31	10
8888591081	Angled cannula	24	38	10
8888591016	Straight cannula	10	13	10
8888591131	Straight cannula	12	14	10
8888591032	Straight cannula	14	15	10
8888591156	Straight cannula	16	16	10
8888591057	Straight cannula	18	17	10
8888591073	Straight cannula	21	22	10
8888591099	Straight cannula	24	25	10

Argyle™ Ferguson Catheter



- Single-use, transparent catheter for ventricle and/or atrium in thin-walled, thermosensitive PVC, latex free

Code	Dimension Fr	Dimensions Outside Diameter	Dimensions Length	Sales Unit
8888590018	18	6mm	51cm	10

Argyle™ Carotid Artery Shunt



- Single-use, transparent shunt in medical grade non-pyrogenic PVC

Code	Length	Sales Unit
8888577775	23cm	10

The code number includes the following shunts: 8Fr (2.7mm); 10Fr (3.3mm); 12Fr (4.0mm); 14Fr (4.7mm)

Cardiothoracic Accessories

Argyle™ Vascular Tourniquet Kit

- Includes 2 sheaths with securing cap, 1 obturator and 2 fixation tapes x 3.2mm x 61cm, latex free


Code	Description	Sales Unit
8888585000	Tourniquet Kit	10



Sentinel Loops™ and Argyle™ Vessel Loops

- Single-use vessel loops in 100% silicon, x-ray opaque, colour coded, latex free

Code	Description	Length	Sales Unit
8888585208	Regular blue	38cm	10
8888585216	Regular red	38cm	10
8888585224	Regular yellow	38cm	10
8888585232	Regular white	38cm	10
8888585240	Small blue	35.5cm	10
8888585257	Small red	35.5cm	10
8888585265	Small yellow	35.5cm	10
8888585273	Small white	35.5cm	10



Penrose Tubing

- Latex tubing
- X-ray opaque
- Indicated for gravity wound drainage

Code	Description I.D.	Description I.D. Length	Sales Unit
8888513002	6mm	30cm	50
8888513200	8mm	30cm	50
8888513408	10mm	30cm	50
8888513606	13mm	30cm	50
8888513804	16mm	30cm	50
8888514000	19mm	30cm	50
8888514208	22mm	30cm	50
8888514406	25mm	30cm	50
8888514604	6mm	45cm	50
8888514802	8mm	45cm	50
8888515007	10mm	45cm	50
8888515205	13mm	45cm	50
8888515403	16mm	45cm	50
8888515601	19mm	45cm	50
8888515809	25mm	45cm	50

