

Silicone Coated Latex Foley Catheters

Curity™ Silicone Coated Latex Foley Catheter



- Silicone coated latex
- Standard/nelaton tip
- Sterile

| Code | Fr/Ch | Length | Sales Unit |
|--|-------|--------|------------|
| 2 Way 3ml Balloon - Paediatric Catheter | | | |
| 2908-02 | 8 | 28cm | 12 |
| 2910-02 | 10 | 28cm | 12 |
| 2 Way 10ml Balloon | | | |
| 1125-02 | 12 | 41cm | 12 |
| 1185-02 | 14 | 41cm | 12 |
| 1213-02 | 16 | 41cm | 12 |
| 1394-02 | 18 | 41cm | 12 |
| 1440-02 | 20 | 41cm | 12 |
| 1575-02 | 22 | 41cm | 12 |
| 1620-02 | 24 | 41cm | 12 |
| 1746-02 | 26 | 41cm | 12 |
| 2 Way 30ml Balloon | | | |
| 26163 | 16 | 41cm | 12 |
| 26183 | 18 | 41cm | 12 |
| 26203 | 20 | 41cm | 12 |
| 26223 | 22 | 41cm | 12 |
| 26243 | 24 | 41cm | 12 |
| 3 Way 30ml Balloon | | | |
| 9718-02 | 18 | 41cm | 12 |
| 9720-02 | 20 | 41cm | 12 |
| 9722-02 | 22 | 41cm | 12 |
| 9724-02 | 24 | 41cm | 12 |

Curity™ Silicone Coated Latex Foley Catheter - Tiemann



- Silicone coated latex
- Tiemann tip
- Sterile

| Code | Fr/Ch | Length | Sales Unit |
|---------------------------|-------|--------|------------|
| 2 Way 10ml Balloon | | | |
| 9812-02 | 12 | 41cm | 12 |
| 9814-02 | 14 | 41cm | 12 |
| 9816-02 | 16 | 41cm | 12 |
| 9818-02 | 18 | 41cm | 12 |
| 9820-02 | 20 | 41cm | 12 |
| 9822-02 | 22 | 41cm | 12 |
| 9824-02 | 24 | 41cm | 12 |