

Alginate Dressings

Curasorb™ Calcium Alginate Dressings



- Absorbs exudate from wound and creates a protective cushioning gel and helps promote an optimal moist healing environment
- Protective gel binds exudate and helps to prevent skin maceration
- Easy, painless dressing removal minimizes tissue trauma and maximizes patient comfort
- Absorbs approximately 20 times dressing's weight
- Sterile
- Plus version available for increased absorbency

Code	Description	Per Box	Ship Case
1 Sterile 1's in Peel-Back Package			
9231	30cm Rope	5	20
9232	5cm x 5cm	10	100
9233	10cm x 10cm	10	50
9236	10cm x 10cm Plus	10	50
9240	10cm x 14cm	10	50
9238	20cm x 10cm	5	50
9239	15cm x 25cm	10	50
9242	30cm x 60cm	5	5
9243	60cm Rope	5	20
9244	91cm Rope	5	20

Curasorb™ ZN Calcium Alginate Dressings



- Absorbent calcium alginate dressing with the added benefit of zinc

Code	Description	Per Box	Ship Case
Sterile 1's in Peel-Back Package			
9351	30cm Rope	5	20
9354	5cm x 5cm	10	100
9355	10cm x 10cm	10	50
9356	10cm x 20cm	5	50

Indications

- Moderate to Heavy Exudating Wounds
- Tunneling Wounds
- Partial and Full-thickness Wounds
- Pressure Ulcers (Stages II, III & IV)