

# Electrodes



ECG Electrodes for Adults  
Reusable ECG Electrodes  
Accessories for Reusable ECG Electrodes  
ECG Electrodes for Children  
Reusable ECG Electrodes for Children  
Neonatal ECG Electrodes  
Kendall™ Infant Disposable ECG Electrodes for Speciality Application  
Multifunction Defibrillation Electrodes  
Kendall DL™ Disposable Lead-Wires  
Kendall DL™ Reuseable Adapters

For general enquiries call +353 14381600  
or fax +353 14073174

[Click here for a full list  
of International contacts](#)

## ECG Electrodes for Adults

## Disposable ECG Electrodes

- With Ag/AgCl sensor and foam backing

Code	Match Code	Description	Quantity Pouch/Carton
Kendall™ Arbo™ Disposable ECG electrode			
31.1946.21	H34SG	Radiolucent, suitable for X-ray, CT and MRI with hydrogel, strong adhesive (72 hour wearable, shower and wash resistant), with integrated abradar for skin preparation, 45mm Ø, with stud adaptor	30/300
31.1925.21	H92SG	With special hydrogel for long-term use, 48 x 34mm, with stud adaptor	50/500
31.1355.21	H135SG	With special hydrogel for long-term use, also suitable for neurology and paediatric applications, 35mm Ø, stud adaptor	50/500
31.1355821	H135TSG	As H135SG, but radiolucent, suitable for CT and MRI	50/500
31.1245.21	H124SG	With special hydrogel for long-term use, also suitable for neurology and paediatric applications, 24mm Ø, stud adaptor	50/500
31.1917.21	H91SSG	As H92SG, but equipped with stress adhesive foam (recommended for long-term ECG monitoring and stress test), 48 x 34mm, with stud adaptor	50/500
31.1915821	H91TSG	As H92SG, but radiolucent, suitable for X-ray, CT and MRI, 48 x 34mm, with stud adaptor	50/500
31.1935.21	H93SG	As H92SG, however only 35 x 26mm	50/500
31.1715.21	H71SG	With integrated cable fixing for motion artefacts reduction, specially suitable for long-term and stress ECG, with stud adaptor	30/300
Kendall™ Medi-Trace™ Disposable ECG Electrode			
31439725	MT Mini 133	30mm Ø, with stud adaptor	3/600
31050522	MT200	36mm Ø, with stud adaptor, packaged on perforated plates	5 x 20/1000
31118733	MT Mini 100	30mm Ø, with stud adaptor, packaged on perforated plates	5 x 20/1000
31013926	MT530	45mm Ø, with stud adaptor	30/600

31.1946.21



31.1925.21



31.1715.21



31050522



## ECG Electrodes for Adults

## Kendall™ Monitoring Disposable ECG Electrodes for Adults

31.1986.21



31.1996.21



- With Ag/AgCl sensor, micropore, cloth, transparent tape or polyester carriers

Code	Match Code	Description	Quantity Pouch/Carton
Kendall™ Arbo™ Disposable ECG Electrode			
31.1986.21	H98SG	Micropore, hydrogel – particularly suitable for sensitive skin, with stud adaptor, 60mm Ø	30/300
31.1996.21	H99SG	Cloth, hydrogel – particularly suitable for sensitive skin, with stud adaptor, 55mm Ø	30/300
Kendall™ Disposable ECG Electrode			
ES40076	H59P	As H49P, however 35 x 22mm	10 x 3ea/450
ES40281	H69P	As H49P, however 22 x 22mm	10 x 3ea/450
Kendall™ Medi-Trace™ ECG Electrodes			
22700	MT700	Kendall™ Medi-Trace™ 700 ECG electrode on perforated plates, each with 10 electrodes with transparent tape backing, 45mm Ø, stud adaptor and hydrogel	5 x 10ea/1000
22733	MT733	Kendall™ Medi-Trace™ 733 ECG electrode, same as Medi-Trace™ 700, but one strip of 3 electrodes	3/600
22735	MT735	Kendall™ Medi-Trace™ 735 ECG electrode, same as Medi-Trace™ 700, but one strip of 5 electrodes	5/600
22750	MT750	Kendall™ Medi-Trace™ 750 ECG electrode, same as Medi-Trace™ 700, but for X-ray, CT and MRI applications	5 x 10ea/1000

31.1345.21



31.1985.21



31.1995.21



- With liquid gel and Ag/AgCl sensor, different backings

Code	Match Code	Description	Quantity Pouch/Carton
31.1345.21	H34LG	Foam backing, strong adhesive (wearable for 72 hours, shower and wash resistant), particularly suitable for stress test and Holter monitoring, 45mm Ø, with stud adaptor	30/300
31.1663.21	H66LG	Foam backing, 55mm Ø, with stud adaptor	30/300
31.1985.21	H98LG	Micropore backing, suitable for sensitive skin, 60mm Ø, with stud adaptor	30/300
31.1995.21	H99LG	Cloth backing, suitable for sensitive skin, 55mm Ø, with popper	30/300

## ECG Electrodes for Adults


## Kendall™ Diagnostic (resting) Disposable ECG Electrodes

- With conductive and adhesive hydrogel and Ag/AgCl sensor

Code	Match Code	Description	Quantity Pouch/Carton
Kendall™ Disposable ECG TAB Electrodes			
31433538	Q-Trace™ 5400	Kendall™ Q-Trace™, full surface adhesive and conductive hydrogel, 23 x 25mm	10 x 10/4000
30807732	Q-Trace™ 5500	Kendall™ Q-Trace™ gold, full surface adhesive and conductive hydrogel, 22 x 31mm	10 x 10/2000
31.1914.21	H914SG	Hydrogel, foam carrier, 48 x 34mm, compatible with 4mm banana plug	30/300
31.1914821	H914TSG	As H914SG, however radiolucent, suitable for X-ray applications	30/300
Adaptor for TAB/SNAP ECG Electrodes			
22660	A34	Universal adaptor for stud and TAB electrodes, adapts to 4mm banana plug	10/-



31433538



3080732



31.1914.21



22660




3390.0400.01

## Reusable ECG Electrodes


## Kendall™ Meditec Chest Suction and Limb Clamp ECG Electrodes

- With Ag/AgCl or silver plated sensor

Code	Description	Quantity Pouch/Carton
2062000330	Kendall™ Meditec Welsh-type chest suction electrode with metal suction cup, 30mm Ø, blue plastic suction ball, with universal adaptor and Ag/AgCl sensor	1
2062000324	Kendall™ Meditec Welsh-type chest suction electrode, as 2062000330, but with 24mm suction cup	1
206700A330	Kendall™ Meditec chest suction electrode with soft plastic suction cup, 30mm Ø, blue plastic suction ball, with universal adaptor, Ag/AgCl sensor	1
206800A330	Kendall™ Meditec chest suction electrode with soft plastic suction cup, 30mm Ø, blue plastic suction ball, with 4mm banana adaptor and Ag/AgCl sensor	1
206800A324	Kendall™ Meditec chest suction electrode with rigid plastic suction cup, 24mm Ø, blue plastic suction ball, with 4mm banana adaptor and Ag/AgCl sensor	1
2060T0030B	Kendall™ Meditec spare suction ball for chest suction electrodes with Ø 30mm cups, blue	1
2011000300	Kendall™ Meditec limb clamp ECG electrodes, IEC colour coded (red, yellow, green, black), with universal adaptor and Ag/AgCl sensor	1 set of 4
2011000307	Kendall™ Meditec limb clamp ECG electrodes, blue coloured, with universal adaptor and Ag/AgCl sensor	1 set of 4
81.101.000	Kendall Arbo™ chest suction ECG electrode with rubber suction cup, 24mm Ø, red rubber ball, with 4mm banana plug adaptor and silver plated sensor	1



206800A324



206700A330