

Reusable ECG Electrodes for Children

Kendall™ Meditec Chest Suction ECG Electrodes and Limb Clamp Electrodes for Children

- With Ag/AgCl sensor

Code	Description	Quantity Pouch/Carton
2062000315	Kendall™ Meditec chest suction electrode Welsh type for children with metal cup, 15mm Ø, with universal adaptor	1



2062000315

Neonatal ECG Electrodes

Disposable ECG Electrodes for Infants

- With conductive and adhesive hydrogel and transparent tape backing Kendall™ Medi-Trace™ Puppy Dog

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
31424780	1041PTS/60cm	AAMI (white, black, green)	1.5mm DIN safety plug	Round	25Ø	Suitable for incubator, radiolucent, repositionable	3/300
31424785	1042PTS/60cm	AAMI (white, black, green)	1.5mm DIN safety plug	Rectangular	25 x 30	Suitable for incubator, radiolucent, repositionable	3/300



31424780



31424785


Kendall™ Medi-Trace™ KittyCat

- With conductive and adhesive hydrogel and cloth backing

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
31424768	1050NPSM/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Round	25Ø	Suitable for incubator, repositionable	3/300
31424743	1051NPSM/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Round	25Ø	Suitable for incubator, repositionable	3/300
31424735	1052NPSM/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Round	30Ø	Suitable for incubator, repositionable	3/300
31424719	1053NPSM/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Round	30Ø	Suitable for incubator, repositionable	3/300




31424768



31424743




31424735



31424719

Neonatal ECG Electrodes

Kendall™ KittyCat ECG Electrodes



30710001



31368072



31086575

- With conductive and adhesive hydrogel and foam backing

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
30710001	4103/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Round	22Ø	Suitable for incubator, repositionable	3/300
31368072	4203/60cm	IEC (red, yellow black)	1.5mm, DIN safety plug	Round	22Ø	Suitable for incubator, repositionable	3/300
31086575	4105/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Oval	20 x 30	Suitable for incubator, repositionable	3/300


Kendall™ Arbo™ Kiddy ECG Electrodes




31.5872.04




31.5273.03



31.5072.03



31.5032.26



31.5033.26

- With conductive and adhesive hydrogel and micropore backing

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
31.5872.04	H87PG/E/60cm	IEC (red, yellow, black)	1.5mm DIN safety plug	Round	30Ø	Suitable for incubator, repositionable	3/120
31.5273.03	H27PG/F/60cm	IEC (red, yellow, black)	1.5mm DIN safety plug	Round	22Ø	Suitable for incubator, repositionable	3/120
31.5072.03	H207PG/F/60cm	IEC (red, yellow, black)	1.5mm DIN safety plug	Round	22 x 30	Suitable for incubator, repositionable	3/120
31.5032.26	H203PG/50cm	White	4mm banana socket	Oval	22 x 30	Suitable for incubator, repositionable	15/120
31.5033.26	H203PGK/10cm	White	4mm banana socket, extra short lead	Oval	22 x 30	Suitable for incubator, repositionable	15/120
31.5072824	H207PGT/90cm	Carbon leads, transparent leads jacketing	1.5mm DIN safety plug	Oval	22 x 30	Suitable for incubator, radiolucent, longer leads, repositionable	3/90
31.5039124	H203PGT/90cm	Carbon leads, transparent leads jacketing	1.5mm DIN safety plug	Oval	22 x 30	Suitable for incubator, radiolucent, longer leads, repositionable	3/90
31.5072.26	H207PGN/150cm	White	1.5mm DIN safety plug	Oval	22 x 30	Extra long lead, suitable for sleep applications, repositionable	10/100

Neonatal ECG Electrodes

Kendall Arbo™ Cutie™ ECG Electrodes

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
31.3876.04	H87H/F/60cm	IEC (red, yellow, black)	1.5mm DIN safety plug	Round	30Ø	Suitable for standard care repositionable	3/120
31.3276.03	H27H/F/60cm	IEC (red, yellow, black)	1.5mm DIN safety plug	Round	22Ø	Suitable for standard care repositionable	3/120
31.3076.26	H207H/60cm	White	1.5mm DIN safety plug	Oval	22 x 30	Suitable for standard care repositionable	15/120



31.3876.04



31.3076.26

Kendall™ Infant Disposable ECG Electrodes for Speciality Application

Kendall™ Medi-Trace™ Limb Band ECG Electrode for Prematures & Speciality Application

- With conductive and adhesive hydrogel and cloth backing


Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
EP30016	LB36-LB/60cm	AAMI (white, black, green)	1.5mm, DIN safety plug	Rectangular	16 x 18	Suitable for incubator, with integrated fastening band, repositionable	3/300




Kendall™ Arbo™ Goldy™ ECG Electrodes for Infants

- With Karaya hydrogel and micropore backing


Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
31.3870.24	H87V/50cm	White	1.5mm DIN safety plug	Round	30Ø	Hypoallergenic, repositionable	15/120
31.3873.25	H87V/F/60cm	IEC (red, yellow, black)	1.5mm DIN safety plug	Round	30Ø	Hypoallergenic, repositionable	3/120
31.4273.26	H27V/60cm	White	1.5mm DIN safety plug	Round	22Ø	Hypoallergenic, repositionable	3/120
31.3835.25	H83V/50cm	White	4mm banana socket	Round	30Ø	Hypoallergenic, repositionable	15/120



31.3870.24



31.3873.25



31.4273.26

Kendall Arbo™ ECG Electrodes for Infants Ungelled

- With micropore backing and NO contact medium

Code	Match Code/ Lead Length	Colour Code	Connector Type	Shape	Size (mm)	Special Features	Qty Sales Pack
31.2870.25	H87L/60cm	White	1.5mm DIN safety plug	Round	30Ø	Ungelled	15/120
31.2830.25	H83L/50cm	White	4mm banana socket	Round	30Ø	Ungelled	15/120