

Perineal Cold Pads

Cool Comfort™ Instant Perineal Cold Pads



- Easy to use and store
- Always ready on shelf
- Non-toxic, non-caustic, non-irritating and latex free
- Cools within a few minutes
- Cooling effect lasts up to 4 hours
- Highly absorbent
- With impermeable backing to prevent fluid leakage
- Self-adhesive tape
- Simple disposal in line with in-house disposal regulations

Code	Description	Sales Unit
5200024	Cool Comfort™ Basic Cold Pad	24
MH20002	Cool Comfort™ with self-adhesive stickers	24

Breast Pads

MotherMates™ Hydrogel Breast Feeding Pads



Features:

- Safe substances used
- Absorbs leaking milk
- Wound healing by creating a physiological wound climate
- Pain relief due to cooling effect
- Can be used for up to 7 days
- Supplied in pairs
- Sterile

Indication:

- To prevent and treat sore and cracked nipples
- Made of polyurethane based hydrophilic hydrogel

Code	Description	Minimum Ordering Qty	Sales Unit
MotherMates™	Hydrogel 7.6cm Ø round		
9403	Packaged in pair	50	50 pairs