

Devon™ Light Handle Systems



The Devon™ Disposable Light Handle System provides an effective way to reduce costs and standardise overhead surgical lights, offering a complete solution to your specific surgical lighting needs. The disposable and sterile system eliminates costly reprocessing and simplifies ordering and inventory management.

Liteglove™ Light Glove



- Flexible light handle cover - convenient and economical

Code	Description	Sales Unit
31140208	Liteglove™ Light handle cover	1 pack/25 box/100 case
31140257	Liteglove™ Light handle cover	2 pack/40 box/160 case
31140216	Liteglove™ Light handle cover	3 pack/48 box/192 case

EZ Handle™ Light Handle



- Rigid, one piece light handle - screws securely in place

Code	Description	Sales Unit
31140133	EZ Handle™ Light handle	1 pack/16 box/128 case
31140125	EZ Handle™ Light handle	2 pack/28 box/112 case
31140117	EZ Handle™ Light handle	3 pack/21 box/84 case

Non-Sterile Disposable Lite Handle - Rigid Handle

Code	Description	Sales Unit
31159455	EZ Handle™ Light handle, non-sterile	300

Litesleeve™ Light Sleeve



- Semi-rigid light handle cover - solid feel

Code	Description	Sales Unit
31140273	Litesleeve™ Light Sleeve	1 pack/70 case
31140240	Litesleeve™ Light Sleeve	2 pack/120 case

Non-Sterile Disposable Lite Sleeve - Semi-Rigid Cover

Code	Description	Sales Unit
31159414	BLite Sleeve, disposable, light handle cover, banded	160

Devon™ Light Handle Systems

Conversion Guide to Light Handles and Light Adaptors

- ULH: Designed for use with the Liteglove™ Light handle cover or Litesleeve™ Light sleeve
- ULA: Designed for use with the EZ Handle™ Light handle



Model Name/ Number	Handle Description	ULH	Reorder Number	ULA	Reorder Number
ALM Surgical					
All ECL and PRC models	Non-centre focus	3600-9	31146866	3660-9	31146569
Angenieux AX4 Major	Centre focus	3600-8	31146585	3650-8	31146452
Angenieux AX14 Major	Centre & non-centre	3600-8	31146585	3650-8	31146452
Daylighter Major Light	Non-centre focus	3600-8	31146585	3650-8	31146452
Daylighter 150, Single Sate.	Non-centre focus	3600-8S	31146858	3650-8S	N/A
Angenieux AX10 Dual	Centre focus	3600-8	31146585	3650-8	31146452
Illuminator III			NA	3650-11	
21st Century series	Centre & non-centre	3600-8	31146585	3650-8	31146452
Surgilite II	Centre & non-centre	3600-8	31146585	3650-8	31146452
All PRX models	Must have ALM install interface EZ handle works with interface	3600-25	31335634		
Prismalix	Same as PRX models	3600-25	31335634		
X-Ten	Same as PRX models	3600-25	31335634		
Medical Illuminations			NA	3650-11	
GCL 500/700		3600-57		3660-57	
Prismatic		3600-9	31146866	3660-9	31146569
751		3600-9	31146866	3660-9	31146569
TRC 7501		3600-9	31146866	3660-9	31146569
Steris					
Gemini	Non-centre focus	3600-2	31146759	3600-2	31146551
Polaris	Non-centre focus	3600-2	31146759	3600-2	31146551
Quantum*	Centre focus	3600-2Q		3660-Q	
Challenge 22	Non-centre focus	3600-2	31146759	3660-2	31146551
SQ240	Non-centre focus	3600-2A		3660-2	31146551
Harmony		3600-38	50001200	3660-38	50001201
L22				3660-2	31146551
Berchtold					
Chromophare 450	Centre focus	3600-22	31146668	3660-12S	31146320
Chromophare 480	Centre focus		NA		NA
Chromophare 570	Centre focus	3600-22	31146668	3660-12S	31146320
Chromophare 950G	Centre focus	3600-95	31146627	3650-95	31146395
Chromophare C-571	Centre focus	3600-12	31146809	3660-12S	31146320
Chromophare D-650	Centre focus	3600-36	50000928	3660-36	50000929
Chromophare	Chromovision	3600-22H	31368627	3660-22H	
Chromophare D-530	Centre focus	3600-36	50000928	3660-36	50000929
Chromophare D-630	Centre focus	3600-36	50000928	3660-36	50000929
Chromophare D-540	Centre focus	3600-36	50000928	3660-36	50000929
Chromophare D-650	Centre focus	3600-36	50000928	3660-36	50000929
Chromophare D-660	Centre focus	3600-36	50000928	3660-36	50000929
Brandon					
Galaxy Ultra		3600-35	50000571		NA
Burton					
11400	Non-centre focus	3600-13	31146700	3660-13	31146429
11410	Non-centre focus	3600-13	31146700	3660-13	31146429
114132	Non-centre focus	3600-13	31146700	3660-13	31146429
Genesis	Centre focus		NA		NA

Devon™ Light Handle Systems

Conversion Guide to Light Handles and Light Adaptors

Model Name/ Number	Handle Description	ULH	Reorder Number	ULA	Reorder Number
Dr. Mach™*					
Quintaflex - R96	Centre focus	3600-14L		NA	
Triaflex	Non-centre focus	3600-14	31146692	3650-14	31146338
Trigenflex - R96	Centre focus	3600-14L		NA	
Uniflex	Centre focus	3600-14	31146692	3650-14	31146338
Getinge Castle™*					
Castle 2000 series		3600-24	31146643	3650-24	
Castle 2023 series		3600-24	31146643	3650-24	
Castle 2300 series	Screw-on	3600-1	31146825	NA	
Castle 2400 series	Screw-on	3600-1	31146825	NA	
Castle 2500-2900 series	Screw-on	3600-3	N/A	3660-3	
Castle 9800 series	Deep-dish	3600-5		3660-5	
Castle Orbiter	Deep-dish	3600-5		3660-5	
Castle Twinstar	Deep-dish	3600-5		3660-5	
Castle Daystar	Deep-dish	3600-5		3660-5	
Castle TV Daystar	Deep-dish	3600-5		3660-5	
Charles Polo Adaptor		3600-5C	31160701	3660-5C	
Sybron 2600		3600-3	N/A	3660-3	
Castle 500/700 series	Snap-on	3600-57		3660-57	
Castle GNM85004	Centre focus snap-on	NA		NA	
Castle TV Daystar	Snap-on	3600-4		3660-4	
Heraeus™*					
Amsterdam	Centre focus	3600-10S	31160081	3660-10S	31160628
Boston	Non-centre focus	3600-10B	31146882	3660-10B	31146270
Frankfurt	Centre focus	3600-10	31146874	3660-10	31146494
Hanau	Centre focus	3600-10H	31146890	3650-10H	31146288
Hanaulux 2002	Centre focus	3600-104	31146619	3650-104	31146403
Hanaulux 2003	Centre focus	3600-104	31146619	3650-104	31146403
Hanaulux 2004	Centre focus	3600-104	31146619	3650-104	31146403
Hanaulux 2005	Centre focus	3600-104	31146619	3650-104	31146403
Hanaulux 2007	Centre focus	3600-107	31146601	3650-107	
Hanaulux 2007 AF	Centre focus	3600-107	31146601	3650-107	
London	Centre focus	3600-10	31146874	3660-10	31146494
Melbourne	Centre focus	3600-10	31146874	3660-10	31146494
Oslo	Centre focus	3600-10	31146874	3660-10	31146494
Satellite-London/Heidleburg	Centre focus	3600-10	31146874	3660-10	31146494
Hanaulux Blue Line	Non-centre focus	3600-104	31146619	3650-104	31146403
Hanaulux 2000		3600-104	31146619	3650-104	31146403
Hannavision		3600-104	31146619	3650-104	31146403
OMI (Midmark)					
Ultralux	Centre focus	3600-23	31146650	3660-23	31181319
Chick 20 series	Non-centre focus	3600-9	31146866	3660-9	31146569
Ultralux M-20	Centre focus	3600-23	31146650	3660-23	31181319
Chick 6100 series			NA		NA
Ultralux 6500 series	Centre focus		NA	3660-23	31181319
Ultralux 6600 series	Non-centre focus use Lite Glove w/OEM			3660-23B	
Ultralux 355 Minor Surgery	Use EZ Handle w/OEM	3600-25	31335634		
Chick II			NA	3660-20	
Ritter 355		3600-25	31335634	NA	

Devon™ Light Handle Systems

Conversion Guide to Light Handles and Light Adaptors

Model Name/ Number	Handle Description	ULH	Reorder Number	ULA	Reorder Number
Stryker™					
Visum			NA	3660-62	50001202
Skytron™					
Nova HA58 & HA58/EL	Green - Major	3600-6S	31146742	3660-6S	31160586
Nova HA75	Green - Major	3600-6S	31146742	3660-6S	31160586
Nova HA75/AR58	Green - Major	3600-6		3660-6	
Nova HA90 & HA90/EL	Green - Major	3600-6		3660-6	
Nova HA100 & HA100/EL	Green - Major	3600-6		3660-6	
58 & 75 & 96	Green - Major	3600-6		3660-6	
5858 & 9656	Green - Major	3600-6S	31146742	3660-6S	31160586
9756 & 9856	Green - Major	3600-6S	31146742	3660-6S	31160586
7558 & 9058 & 10058	Green - Major	3600-6		3660-6	
Nova HA58/EL	Green - Minor	3600-6S	31146742	3660-6S	31160586
Nova HA75/EL	Green - Minor	3600-6		3660-6	
Nova HA75/AR58	Green - Minor	3600-6S	31146742	3660-6S	31160586
Nova HA76/AR58	Green - Minor	3600-6S	31146742	3660-6S	31160586
Nova HA90 & HA100/EL	Green - Minor	3600-6S	31146742	3660-6S	31160586
5858 & 7558	Green - Minor	3600-6S	31146742	3660-6S	31160586
9058 & 10058	Green - Minor	3600-6S	31146742	3660-6S	31160586
9656 & 9756 & 9856	Green - Minor	3600-7	31146833	3660-7	
Nova HA58 & HA58/EL	Tan - Major	3600-7N	31146841	3660-7N	31146460
Nova HA75 & HA 90	Tan - Major	3600-7N	31146841	3660-7N	31146460
Nova HA90 & HA100/EL	Tan - Major	3600-7N	31146841	3660-7N	31146460
Nova HA100	Tan - Major	3600-7N	31146841	3660-7N	31146460
Nova HA75/AR58	Tan - Major	3600-6		3660-6	
Infinity IN30 & IN22 series	Tan - Major	3600-7N	31146841	3660-7N	31146460
58 & 75	Tan - Major	3600-7N	31146841	3660-7N	31146460
5858 & 7558	Tan - Major	3600-7N	31146841	3660-7N	31146460
9058 & 10058	Tan - Major	3600-7N	31146841	3660-7N	31146460
96 & 9656	Tan - Major	3600-7	31146833	3660-7	
9756 & 9856	Tan - Major	3600-7	31146833	3660-7	
Nova HA58/EL	Tan - Minor	3600-7N	31146841	3660-7N	31146460
Nova HA75/EL	Tan - Minor	3600-7N	31146841	3660-7N	31146460
Nova HA76/AR58	Tan - Minor	3600-7N	31146841	3660-7N	31146460
Nova HA75/AR58	Tan - Minor	3600-7N	31146841	3660-7N	31146460
Nova HA90 & HA100/EL	Tan - Minor	3600-7N	31146841	3660-7N	31146460
Infinity IN30 & IN22 series	Tan - Minor	3600-7N	31146841	3660-7N	31146460
5858 & 7558	Tan - Minor	3600-7N	31146841	3660-7N	31146460
9058 & 10058	Tan - Minor	3600-7N	31146841	3660-7N	31146460
9756 & 9656 & 9856	Tan - Minor	3600-7	31146833	3660-7	
5656	Both	3600-7	31146833	3660-7	
Excel KR56	Both	3600-7	31146833	3660-7	
Excel SH60	Both	3600-7	31146833	3660-7	
Excel SH61	Both	3600-7	31146833	3660-7	
Excel SH63	Both	3600-7	31146833	3660-7	
Universe II		3600-72	31146635	3660-72	
Skylux		3600-17S	31430195	3660-17S	
Universe 60/63		3600-7	31146833	3660-7	
Infinity IF series		3600-7E	31146726	3660-7E	31146478
Infinity 5400 series		3600-7E	31146726	3660-7E	31146478
Infinity IF series	With Camera	3600-7H		3660-7H	
Universe w/centre focus			NA		NA
Skytron 16			NA		NA
Stellar Series		3600-7E	31146726	3660-7E	31146478
Hill-Rom™					
Prima	All light controls are in the handle thus no adaptor needed use the LiteGlove				