

Devon™ Mini Kits

Mini Kits

Code	Description	Sales Unit
31141784	K-1960-S Mini-Kit	64
	Litegloves™ (31140257)	2
	Foam block needle counter (31142352)	1
	Tipolisher (31142790)	1
	Suture bag (31142691)	1
	Surgical marker, ruler and label (31145777)	1



OR Essentials

Cautery Tipolisher

Safely and effectively cleans electro-surgical pencil tips, strong adhesive backing and X-ray detectable.

Code	Description	Sales Unit
31142790	Cautery Tipolisher (5 x 5cm)	25 box/100 case
31142717	Cautery Tipolisher (5 x 5cm)	80 box/240 case



Micro Tip Wipe

A sterile, lint and particulate-free foam wipe for cleaning the fine tips of micro instruments and lenses.

Code	Description	Sales Unit
31142535	Micro Tip Wipe	25 box/100 case

Dev-O-Loops Surgical Vessel Loops

Made of silicone; radiopaque bands that find practical application in any situation where arteries, veins, nerves or tendons must be looped for identification, retraction or occlusion. Available in four colours (red, white, yellow, blue) and two sizes.

Code	Description	Colour	Size	Sales Unit
31145728	Dev-O-Loops	White	Mini	20 box/200 case
31145710	Dev-O-Loops	Red	Mini	20 box/200 case
31145702	Dev-O-Loops	Blue	Mini	20 box/200 case
31145660	Dev-O-Loops	Yellow	Mini	20 box/200 case
31145686	Dev-O-Loops	White	Maxi	20 box/200 case
31145744	Dev-O-Loops	Red	Maxi	20 box/200 case
31145678	Dev-O-Loops	Blue	Maxi	20 box/200 case
31145694	Dev-O-Loops	Yellow	Maxi	20 box/200 case




Dev-O-Clamps Drape Clamps

Available as:

- Metal reusable
- Sterile disposable
- Non-sterile disposable

Code	Description	Sales Unit
31140315	Disposable clamp, non-sterile	Bulk/50 case
31140299	Disposable clamp, non-sterile	Bulk/100 case
31140281	Disposable clamp, sterile	1 pack/10 box/80 case
31140307	Disposable clamp, sterile	2 pack/30 box/120 case



Tube Holder



- Organises and securely holds tubes and cords in place
- Accommodates any size tube
- Available with removable or fixed touch fasten strap
- Sterile, single-use

Code	Description	Sales Unit
31142667	Tube holder with removable touch fasten strap	25 box/100 case
31142659	Tube holder with removable touch fasten strap	12 box/96 case
31142675	Tube holder with fixed touch fasten strap	25 box/100 case

OR Utility Products

Code	Description	Sales Unit
31160297	Tube holder, bulk, non-sterile	200
31142675	Tube holder with fixed touch fasten strap (25 per box/4 boxes per case)	100
31160370	Tipolisher, large 5 x 5cm, bulk, sterile	600
33160099	Micro Tip Wipe, bulk, sterile	300
31160131	Scope adjustment bag, bulk, non-sterile	1000
31350757	Operating room suture bag, bulk, sterile	2000
31160222	Blade Shield - scalpel holder	250

Dev-O-Pouch Instrument Pouch



Stores a variety of instruments connected to tubes and cords.

Code	Description	Sales Unit
31142550	Instrument pouch (double pocket)	20 box/80 case

Suture Bag



Adheres securely to any surface for convenient, quick access.

Code	Description	Sales Unit
31142691	Suture bag	25 box/100 case

Scope Adjustment Bag



Sterile microscope cover allows surgeon to adjust microscope during procedure.

Code	Description	Sales Unit
31142543	Scope adjustment bag	120 box/480 case