

Surgical Suction



Argyle™ Surgical Suction
Cannulas

Vital Vue™ Surgical
Instrument

Tracheobronchial Suction
Catheters

Sep-T-Vac™ Collection System
Connecting Tubing

For general enquiries call +353 14381600
or fax +353 14073174

[Click here for a full list
of International contacts](#)

Argyle™ Surgical Suction Cannulas

A comprehensive range of suction cannulas designed to meet every need in the Operating Room.

Argyle™ Yankauers



- Complete range of rigid and flexible single use yankauers for suction of fluids



Argyle™ Flexible Yankauers

- Sterile, single use
- Flexible yankauers can be shaped by the user to modify angle and adapt to specific operative sites

Code	Description	Tip Trol	Size	Length	Sales Unit (Box/packaging)
Argyle™ Flexible Yankauer for Paediatric or Dental Applications					
1180113001	No side holes	No	16Fr	16cm	100
1180113050	4 side holes	Yes	16Fr	16cm	100
1180113100	No side holes	Yes	16Fr	16cm	100
Argyle™ Flexible Yankauer					
8888504001	Fine tip	No	12Fr	25cm	50
8888504019	Fine tip	Yes	12Fr	25cm	50
8888501007	Medium tip	No	18Fr	25cm	50
8888501023	Medium tip	Yes	18Fr	25cm	50
8888501015	Medium tip with sump	No	18Fr	25cm	50
8888501031	Medium tip with sump	Yes	18Fr	25cm	50
1180501205	Medium tip with Sure-Grip	No	18Fr	25cm	50
8888502005	Large tip	No	22Fr	25cm	50
8888502013	Large tip with sump	Yes	22Fr	25cm	50
1180502039	Large tip with Sure-Grip	No	22Fr	25cm	50
1180504018	Flat tip	No	22Fr	25cm	50
1180503003	Extra large tip	No	30Fr	25cm	50
Argyle™ Flexible Yankauer Sets with Non-Conductive Flexible Tubing					
1180504613	Fine tip 5mm dia tube	No	12Fr	2.2m**	25
1180504639	Fine tip 5mm dia tube	No	12Fr	3.0m**	25
1180504621	Fine tip 5mm dia tube	No	12Fr	4.0m**	25
1180501619	Medium tip 7mm dia tube	No	18Fr	2.2m**	25
1180501635	Medium tip 7mm dia tube	No	18Fr	3.0m**	25
1180501700	Medium tip 7mm dia tube	Yes	18Fr	3.0m**	25
1180502617	Large tip 7mm dia tube	No	22Fr	2.2m**	25
1180502625	Large tip 7mm dia tube	No	22Fr	3.0m**	25
1180502618	Large tip 7mm dia tube	No	22Fr	2.2m**	25
1180502658	Large tip 7mm dia tube	No	22Fr	3.0m**	25



Argyle™ Surgical Suction Cannulas

Argyle™ Flexible Yankauers (continued)

Code	Description	Tip Trol	Size	Length	Sales Unit (Box/packaging)
Argyle™ Flexible Yankauer Sets with Conductive Flexible Tubing					
8888501502	Medium tip 7mm dia tube	No	18Fr	1.8m**	50
8888501528	Medium tip with sump 7mm dia tube	Yes	18Fr	1.8m**	50
8888502518	Large tip with sump 7mm dia tube	Yes	22Fr	1.8m**	50

Argyle™ Rigid Yankauers



- Sterile, single use
- Shatter resistant PVC enables instrument to be used for moderate organ retraction or tissue dissection

Code	Description	Tip Trol	Size	Length	Sales Unit (Box/packaging)
8888505115	Large tip	No	22Fr	25cm	50
8888505123	Large tip	Yes	22Fr	25cm	50
8888505016	Bulb tip	No	34Fr	26cm	50
8888505024	Bulb tip	Yes	34Fr	26cm	50

Argyle™ Selec-Trol™ Yankauers



- Sterile, single use
- Available with large open tip, bulb tip or fine tip
- Open and bulb tips are made from transparent high quality PVC with lateral eyes, rounded distal ends and bevelled edges for controlled atraumatic suction, even on the most delicate of tissues

Code	Description	Tip Trol	Size	Length	Sales Unit (Box/packaging)
1180505172	Large tip	On-off switch	24Fr	27cm	20
1180505164	Bulb tip	On-off switch	32Fr	27cm	20

Argyle™ Poole Suction Instruments



- Sterile, single use
- Two-piece fixed construction: transparent inner and outer lumen with fingertip suction control

Code	Description	Tip Trol	Size	Length	Sales Unit Box/packaging
Argyle™ Poole Instruments					
8888508002	Flexible Poole	-	34Fr	27cm	50
8888509703	Rigid Poole	-	14-28Fr	25cm	50
Argyle™ Poole Sets					
1180508515	Flexible Poole with 5mm dia non-cond tube	-	34Fr	2.2m**	25
8888509729	Rigid Poole with 5mm dia cond tube	-	14-28Fr	1.8m**	50

Vital Vue™ Surgical Instrument

3-in-1 Design

The 3-in-1 design of the Vital Vue™ Surgical instrument provides the surgeon with convenient and immediate access between illumination irrigation and aspiration.

Tip Options

The Vital Vue™ Surgical instrument is available in five distinct distal procedural tips, allowing the instrument to be used in a number of surgical procedures.

Vital Vue™ Surgical instrument - Bulb Tip

Used in many gynaecological and plastic surgery cases as well as other shallow procedures where a bulb tip is appropriate.



Vital Vue™ Surgical instrument - Extended Bulb Tip

For lower abdominal and pelvic cases where extra reach is required and a bulb tip is appropriate.



Vital Vue™ Surgical instrument - Extended Slender Tip

The tip of choice for surgeons who prefer precise suction with a slender, 1.35cm profile tip. Especially applicable in cases which have a vascular, ductal or tubular nature.



Vital Vue™ Surgical instrument - Extended Metal Tip

The 12CH, 5cm extended metal tip is ideal where reach and precise aspiration is required, such as oral surgery and some vascular cases.




Vital Vue™ Surgical instrument - Extra Long Extended Metal Tip

Designed specifically for orthopaedic cases, this 12CH, 8cm metal tip is ideal where extra extension is necessary particularly in hip arthroplasty.



Code	Description	Sales Unit (Box/packaging)
Vital Vue™ Surgical Instrument		
8886827906	Vital Vue™ Surgical instrument bulb tip	6
8886828006	Vital Vue™ Surgical instrument extended bulb tip	6
8886828106	Vital Vue™ Surgical instrument extended slender tip	6
8886828206	Vital Vue™ Surgical instrument extended metal tip	6
8886828306	Vital Vue™ Surgical instrument extra long extended metal tip	6
Vital Vue™ Light Power Source		
8886828421	GB Vital Vue™ Surgical instrument light power source X1 (3-pin UK plug)	1
8886828422	ES Vital Vue™ Surgical instrument light power source X1 (2-pin EU plug)	1
8886828423	DE Vital Vue™ Surgical instrument light power source X1 (2-pin EU plug)	1
8886828424	IT Vital Vue™ Surgical instrument light power source X1 (2-pin EU plug)	1
8886828425	FR Vital Vue™ Surgical instrument light power source X1 (2-pin EU plug)	1
F090481	European 2-pin plug for 220v English light power source	1



- A. Light plug
- B. Suction connector
- C. Spike
- D. Irrigation button
- E. Suction bleed hole
- F. Irrigation button retainer clip
- G. Bright obturator (stylet)