

Tandem-Cath™ Dual Lumen Catheter Repair Kits



| Code | Description | Sales Unit |
|------------|---|------------|
| 8888200000 | Tandem-Cath™ repair kit without sideholes | 1 |
| 8888200006 | Tandem-Cath™ repair kit with 6 sideholes | 1 |

Swan Neck Tenckhoff Peritoneal Catheters



| Code | Description | Sales Unit |
|-------------|--|------------|
| 8888412007 | Swan Neck Tenckhoff, 2 cuff left, 43cm | 1 |
| 8888412015S | Swan Neck Tenckhoff, 2 cuff right, 43cm | 1 |
| 8888412023 | Swan Neck Tenckhoff 2 cuff long left, 55.5cm | 1 |

Swan Neck Curl Cath Peritoneal Catheters



| Code | Description | Sales Unit |
|------------|--|------------|
| 8888413807 | Swan Neck Curl Cath 2 cuff left, 62.5cm | 1 |
| 8888413815 | Swan Neck Curl Cath 2 cuff right, 62.5cm | 1 |

Swan Neck Presternal Peritoneal Catheters



| Code | Description | Sales Unit |
|------------|--|------------|
| 8888414011 | Swan Neck Presternal Catheter 2 cuff, 112.8cm (Swan Neck: 60.3cm + 52.5 curl) | 1 |

Swan Neck Missouri Peritoneal Catheters



| Code | Description | Sales Unit |
|------------|---|------------|
| 8888412601 | Swan Neck Missouri 2 cuff left, 44.5cm | 1 |
| 8888412619 | Swan Neck Missouri 2 cuff right, 44.5cm | 1 |

Swan Neck Curl Cath Missouri Peritoneal Catheters

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888413401 | Swan Neck Curl Cath Missouri 2 cuff left, 62cm | 1 |
| 8888413419 | Swan Neck Curl Cath Missouri 2 cuff right, 62cm | 1 |



Swan Neck Oreopoulos-Zellerman (Toronto Western) Peritoneal Catheters

| Code | Description | Sales Unit |
|------------|------------------------------------|------------|
| 8888413005 | Swan Neck O-Z 2 cuff left, 44.5cm | 1 |
| 8888413013 | Swan Neck O-Z 2 cuff right, 44.5cm | 1 |



Oreopoulos-Zellerman (Toronto Western) Missouri Tenckhoff Peritoneal Catheters

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888411009 | Oreopoulos-Zellerman 1 bead, 2 disk, 2 cuff, 41cm | 1 |



Oreopoulos-Zellerman (Toronto Western) Tenckhoff Peritoneal Catheters

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888411025 | Oreopoulos-Zellerman 2 disk, 2 cuff, 41cm | 1 |



Quinton™ Curl Cath Peritoneal Catheters - Two Cuffs



- Beta-Cap™ adapter, cap and clamp
- Curled tubing with numerous inflow/outflow holes
- Can be inserted percutaneously or surgically
- Made of translucent, medical-grade silicone rubber
- Radiopaque stripe
- Can be extended or repaired if necessary
- Sterile

| Code | Description | Sales Unit |
|------------|-------------|------------|
| 8811313015 | Curled 57cm | 1 |
| 8811313010 | Curled 62cm | 1 |
| 8888411702 | Curled 62cm | 1 |

Quinton™ Curl Cath Peritoneal Catheters - One Cuff



| Code | Description | Sales Unit |
|------------|-------------|------------|
| 8811313014 | Curled 57cm | 1 |
| 8811313013 | Curled 60cm | 1 |
| 8888411710 | Curled 62cm | 1 |

Quinton™ Curl Cath Peritoneal Catheters - Loose Cuffs



| Code | Description | Sales Unit |
|------------|-------------|------------|
| 8811313009 | Curled 62cm | 1 |

Quinton™ Tenckhoff Peritoneal Catheters - Two Cuffs



- Beta-Cap™ adapter, cap and clamp
- Can be trimmed to fit most patient sizes
- Sterile

| Code | Description | Sales Unit |
|------------|------------------------------------|------------|
| 8888411405 | Tenckhoff, 2 cuff, 41cm | 1 |
| 8810888003 | Tenckhoff, 2 cuff, 42cm | 1 |
| 8810888012 | Tenckhoff, universal, 2 cuff, 47cm | 1 |
| 8888411447 | Tenckhoff, 2 cuff, 47cm | 1 |

Quinton™ Tenckhoff Peritoneal Catheters - One Cuff

| Code | Description | Sales Unit |
|------------|--|------------|
| 8888411439 | Tenckhoff, 1 cuff, 41cm | 1 |
| 8814843001 | Tenckhoff, 1 preperitoneal cuff, 41cm | 1 |
| 8814843002 | Tenckhoff, universal, 1 preperitoneal cuff, 46cm | 1 |
| 8810889003 | Tenckhoff, 1 subcutaneous cuff, 42cm | 1 |
| 8810889011 | Universal, 1 subcutaneous cuff, 46cm | 1 |
| 8888411421 | Tenckhoff, 1 cuff, 47cm | 1 |



Quinton™ Tenckhoff Peritoneal Catheters - Loose Cuffs

| Code | Description | Sales Unit |
|------------|-----------------|------------|
| 8810887003 | Tenckhoff, 42cm | 1 |



Quinton™ Pediatric/Infant Peritoneal Catheters

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888413100 | Infant Swan Neck Curl Cath, 2 cuff, 38.9cm | 1 |
| 8888413101 | Pediatric Swan Neck Curl Cath, 2 cuff, 43cm | 1 |
| 8888413102 | Adolescent Swan Neck, Curl Cath, 2 cuff, 59cm | 1 |
| 8888414813 | Pediatric Swan Neck Curl Cath, 2 cuff, left, 42cm | 1 |
| 8888415307 | Pediatric Swan Neck Curl Cath, Missouri, 2 cuff, left, 42.5cm | 1 |
| 8888410506 | Pediatric Swan Neck Tenckhoff, 2 cuff, left, 37.5cm | 1 |
| 8888410514 | Pediatric Swan Neck Tenckhoff, 2 cuff, right, 37.5cm | 1 |
| 8815677001 | Pediatric Tenckhoff Curl Cath, 1 cuff, 39cm | 1 |
| 8888414508 | Pediatric Tenckhoff Curl Cath, 1 cuff, 39.25cm | 1 |
| 8888414201 | Pediatric Tenckhoff, 2 cuff, 30cm | 1 |
| 8812329002 | Pediatric Tenckhoff, 2 cuff, 31cm | 1 |
| 8888414227 | Pediatric Tenckhoff, 2 cuff, 31cm | 1 |
| 8810890014 | Pediatric Tenckhoff, 2 cuff, 32cm | 1 |
| 8810890003 | Pediatric Tenckhoff, 2 cuff, 37cm | 1 |
| 8888414219 | Pediatric Tenckhoff, 1 cuff, 30cm | 1 |
| 8812329001 | Pediatric Tenckhoff, 1 subcutaneous cuff, 31cm | 1 |
| 8888414235 | Pediatric Tenckhoff, 1 cuff, 31cm | 1 |
| 8812321001 | Pediatric Tenckhoff, 1 subcutaneous cuff, 37cm | 1 |



Swan Neck Curl Cath Kits



| Code | Description | Sales Unit |
|------------|--|------------|
| 8888413823 | Swan Neck Curl Cath kit, 2 cuff, left, 62.5cm | 1 |
| 8888413831 | Swan Neck Curl Cath kit, 2 cuff, right, 62.5cm | 1 |

Quinton™* Curl Cath Catheter Kits



| Code | Description | Sales Unit |
|------------|---|------------|
| 8817278007 | Quinton™* Curl Cath, 2 cuff, 57cm | 5 |
| 8817278006 | Quinton™* Curl Cath, 2 cuff, 62cm | 5 |
| 8817278010 | Quinton™* Curl Cath, 1 preperitoneal cuff, 60cm | 5 |
| 8817278008 | Quinton™* Tenckhoff, 2 cuff, 42cm | 5 |
| 8817278001 | Quinton™* Tenckhoff Universal, 1 preperitoneal cuff, 46cm | 5 |

Quinton™* Peri-Patch™ Repair Kit



| Code | Description | Sales Unit |
|------------|------------------------|------------|
| 8888285001 | Peri-Patch™ repair kit | 1 |

Quinton™* Single Lumen Right Atrial Catheters

- Made of soft, biocompatible silicone rubber
- Radiopaque
- Can be trimmed to appropriately position tip in right atrium
- Reinforced tubing to protect the catheter when clamping
- Felt cuff to anchor catheter subcutaneously
- Extensions imprinted with name and priming volume
- Can be inserted percutaneously using a Pull-Apart sheath/dilator or surgically implanted
- In-line clamp and luer-lock adapter
- Atraumatic tunneling stylet
- Sterile



| Code | Overall Length | O.D. mm | I.D. mm | Initial Priming Volume (per lumen) | Correction Factor cc/cm | Ship Case |
|------------|----------------|---------|---------|------------------------------------|-------------------------|-----------|
| 8814518001 | 94cm | 2.2 | 1.0 | 1.0cc | 0.0081 | 1 |

Quinton™* Raaf Dual Lumen Right Atrial Catheter Repair Kits

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888890001 | 2.2mm Raaf dual lumen right atrial catheter | 1 |



Quinton™* Single Lumen Right Atrial Catheter Repair Kits

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888665001 | 2.2mm Single lumen right atrial catheter repair kit | 1 |



Quinton™ Raaf Dual Lumen Right Atrial Catheters



- Made of soft, biocompatible silicone rubber
- Radiopaque
- Clear extensions
- Can be trimmed to appropriately position tip in right atrium
- Round lumina in 2.2 and 3.2mm catheters
- Large D-shaped lumina in 4.5mm catheter
- Can be inserted percutaneously using a Pull-Apart sheath/dilator or surgically implanted
- Felt cuff anchors the catheter subcutaneously
- Extensions imprinted with name and priming volume
- Two injection sealing caps
- In-line clamps and luer-lock adapters
- Atraumatic tunneling stylet
- Sterile

| Code | Catheter | Overall Length | O.D. mm | I.D. mm | Initial Priming Volume (per lumen) | Correction Factor cc/cm | Ship Case |
|------------|----------|----------------|-----------|-----------|------------------------------------|-------------------------|-----------|
| 8814850001 | 2.2mm | 90cm | 2.2 x 3.7 | 1.1 | 1.2cc | 0.0089 | 1 |
| 8815331001 | 3.2mm | 102cm | 2.8 x 4.9 | 1.5 | 2.2cc | 0.0188 | 1 |
| 8814525001 | 4.5mm | 95cm | 4.5 | 3.0 x 1.2 | 2.5cc | 0.0233 | 1 |

Quinton™ 3.9mm Triple Lumen Right Atrial Catheter



- Made of soft, biocompatible silicone rubber
- Clear extensions
- Can be trimmed to appropriately position tip in right atrium
- Larger 1.2mm lumen has red adapter
- Felt cuff anchors catheter subcutaneously
- Can be inserted percutaneously through a 12Fr Pull-Apart sheath/dilator or surgically implanted
- Three injection sealing caps
- In-line clamps and color-coded, luer-lock adapters
- Sterile

| Code | Lumen | Initial Priming Volume | Correction Factor cc/cm | Sales Unit |
|------------|---------------|------------------------|-------------------------|------------|
| 8815740001 | White adapter | 1.0cc | 0.0081 | 1 |
| | Blue adapter | 1.0cc | 0.0081 | |
| | Red adapter | 1.25cc | 0.0138 | |

Quinton™ 3.9mm Triple Lumen Right Atrial Catheter Repair Kit



| Code | Description | Sales Unit |
|------------|--|------------|
| 8888824001 | 3.9mm Triple lumen, right atrial catheter repair kit | 1 |

Quinton™ 3.9mm Triple Lumen Right Atrial Catheter

- Made of soft, biocompatible silicone rubber
- Clear extensions
- Can be trimmed to appropriately position tip in right atrium
- Larger 1.2mm lumen has red adapter
- Felt cuff anchors catheter subcutaneously
- Can be inserted percutaneously through a 12Fr Pull-Apart sheath/dilator or surgically implanted
- Three injection sealing caps
- In-line clamps and color-coded, luer-lock adapters
- Sterile



| Code | Lumen | Initial Priming Volume | Correction Factor cc/cm | Sales Unit |
|------------|---------------|------------------------|-------------------------|------------|
| 8815740001 | White adapter | 1.0cc | 0.0081 | 1 |
| | Blue adapter | 1.0cc | 0.0081 | |
| | Red adapter | 1.25cc | 0.0138 | |

Quinton™ 3.9mm Triple Lumen Right Atrial Catheter Repair Kit

| Code | Description | Sales Unit |
|------------|--|------------|
| 8888824001 | 3.9mm Triple lumen, right atrial catheter repair kit | 1 |



Quinton™ Pull-Apart Introducer Sets

| Code | Description | Sales Unit |
|------------|-------------------------------|------------|
| 8815544004 | 7Fr Pull-Apart introducer set | 5 |
| 8815544003 | 7Fr Pull-Apart introducer set | 5 |
| 8815544002 | 7Fr Pull-Apart introducer set | 5 |
| 8815544001 | 7Fr Pull-Apart introducer set | 5 |
| 8815544010 | 7Fr Pull-Apart introducer set | 5 |
| 8815544013 | 7Fr Pull-Apart introducer set | 5 |



8Fr / 10Fr Dilator



| Code | Description | Sales Unit |
|------------|---------------------------|------------|
| 8815072001 | 8Fr Dilator, length 15cm | 5 |
| 8815642001 | 10Fr Dilator, length 15cm | 5 |

J / Straight Stainless Steel Guidewire



- Guides dilator or catheter into a vein or peritoneal cavity
- 2 soft ends, 1 in a J configuration
- Tip straightener attached

| Code | Description | Sales Unit |
|------------|-----------------------------------|------------|
| 8813796001 | J/Straight 0.038" guidewire, 70cm | 10 |
| 8817231001 | J/Straight 0.035" guidewire, 70cm | 10 |

Quinton™* Y-Adapter

- Disposable
- Converts single extension to dual
- Sterile

| Code | Description | Sales Unit |
|------------|----------------------|------------|
| 8818782001 | Disposable Y-Adapter | 25 |



Quinton™* - Injection Sealing Caps

- Luer lock fitting with a natural rubber stopper
- Provides a self-sealing injection site for administering medication or drawing blood
- Universal adapter
- May be used with both haemodialysis and peritoneal dialysis catheters
- Sterile

| Code | Description | Sales Unit |
|------------|------------------------|------------|
| 8813791001 | Injection sealing caps | 50 |



Quinton™* - Peritoneal Tunneling Stylet

- Made of pliable stainless steel
- Can be bent to create customised, curved tunnel
- Sterile

| Code | Description | Sales Unit |
|------------|------------------------|------------|
| 8817791001 | 5Fr x 23cm (9") stylet | 5 |



Quinton™* - 2 Part Titanium Connector

| Code | Description | Sales Unit |
|------------|-------------------------------|------------|
| 8888415604 | 2 Part titanium connector | 1 |
| 8888425769 | Paediatric titanium connector | 1 |



Quinton™* - Beta-Cap™ Adapter

- Made of ULTEM™ resin
- For use with all Quinton™ peritoneal catheters
- Sterile

| Code | Description | Sales Unit |
|------------|-------------------|------------|
| 8814661001 | Beta-Cap™ adapter | 25 |



Quinton™ Introducer Needle



- Stainless steel
- Channels guidewire into a vein or peritoneal cavity
- Accommodates guidewires up to 0.038" diameter
- Luer lock hub
- Sterile

| Code | Description | Sales Unit |
|------------|-----------------------------|------------|
| 8813795001 | 18G x 7cm Introducer needle | 10 |

Quinton™ Faller Tunneling Stylet



| Code | Description | Sales Unit |
|------------|-------------------------|------------|
| 8888415679 | Faller tunneling stylet | 1 |

Quinton™ Peritoneal Catheter Straightening Stylet



- Stainless steel stylet
- Used to keep the catheter straight and manoeuvrable during insertion into the peritoneum
- Autoclavable
- Nonsterile

| Code | Description | Sales Unit |
|------------|--|------------|
| 8888415646 | Insertion/straightening stylet, 44.5cm | 1 |
| 8888415661 | Insertion/straightening stylet, 61cm | 1 |

Quinton™ Titanium Catheter Extender



| Code | Description | Sales Unit |
|------------|----------------------------------|------------|
| 8888415612 | Adult titanium catheter extender | 1 |

Quinton™ Beta-Cap™ Clamp

- For use with all Quinton™ peritoneal catheters
- Sterile

| Code | Description | Sales Unit |
|------------|-----------------|------------|
| 8810805001 | Beta-Cap™ clamp | 25 |



Medical Grade Silicone Adhesive Type A

- Used to bond felt cuffs to catheter
- Also for use with Quinton™ repair kits to extend or repair catheters
- Sterile

| Code | Description | Sales Unit |
|------------|---------------------------|------------|
| 8810807001 | Silicone adhesive, type A | 1 |



Felt Cuffs

- Anchor catheter and help establish a natural infection barrier
- Sterile

| Code | Description | Sales Unit |
|------------|----------------------|------------|
| 8811308001 | Felt cuffs, 1cm long | 10 |



Quinton™ Catheter Stencil Swan Neck - Tenckhoff

| Code | Description | Sales Unit |
|------------|---|------------|
| 8888412304 | Swan Neck Tenckhoff stencil, 2 cuff, left/right | 1 |



Quinton™ Catheter Stencil Swan Neck - Missouri

| Code | Description | Sales Unit |
|------------|--|------------|
| 8888412312 | Swan Neck Missouri stencil, 2 cuff, left/right | 1 |

